




Single-Use Ophthalmology Instruments




Contents

Forceps	001-012
Holder	013-016
Scissors	017-025
Probes, hook & spatulas	025-042
Speculums	043-050
Markers	050-051
Vitreoretinal Instruments	052-054
Cannulas	055-058
Instruments Sets	059-060


FORCEPS




SF-11111-2 Toothed Forceps
85mm, 0.12mm 1 x 2 teeth, straight with tying platforms, flat handle




SF-11112-2 Toothed Forceps
105mm, 0.12mm 1 x 2 teeth, straight with tying platforms, flat handle




SF-11141-2 Toothed Forceps
105mm, 0.12mm 1 x 2 teeth, straight with tying platforms, round handle



SF-11131-2 Toothed Forceps
105mm, 0.12mm 1 x 2 teeth, straight with tying platforms, flat handle




SF-11117-2 Colibri Conjunctiva Toothed Forceps : 85mm, 0.12mm 1 x 2 teeth, curved with tying platforms, flat handle




Colibri Toothed Forceps : 0.12mm 1 x 2 teeth, curved with tying platforms, flat handle


SF-11113-2 85mm
SF-11123-2 110mm




SF-11114-2 Mcpherson Toothed Forceps
85mm, 0.12mm 1 x 2 teeth, angled with tying platforms



SF-11159-1 Adsons Forceps: 125mm, 0.3mm, 1 x 2 teeth, straight without platforms




SF-11061-1 Hoskins / Notched Forceps: 85mm, 0.25mm notched, straight with tying platforms,




SF-11061-2 Hoskins / Notched Forceps
105mm, 0.25mm notched, straight with tying platforms,

FORCEPS




SF-11081-1 Hoskins / Notched Forceps
105mm, 0.25mm notched, straight with tying platforms,



SF-11084-2 Hoskins / Notched Forceps:
85mm, 0.25mm notched, curved with tying platforms,




SF-11101-1 Pierse-Colibri Corneal Forceps
85mm, 0.25mm notched, curved with tying platforms,



SF-11174-1 Capsulorhexis Forceps
85mm, Curved shafts, very delicate grasping tips


SF-11174-1 1.8 MICS Capsulorhexis Forceps
1.8mm incision



SF-11174-2 Capsulorhexis Forceps
105mm, Curved shafts, very delicate grasping tips

SF-11174-2 1.8 MICS Capsulorhexis Forceps
1.8mm incision


FORCEPS



SF-11223-1 Capsulorhexis Forceps
105mm, Curved shafts, very delicate grasping tips




SF-11223-1 1.8 MICS Capsulorhexis Forceps
1.8mm incision



SF-11201C-1 Capsulorhexis Forceps
105mm, Curved shafts, very delicate grasping tips




SF-11201C-1 1.8 MICS Capsulorhexis Forceps
1.8mm incision



SF-11226-1 Capsulorhexis Forceps
120mm, Curved shafts, very delicate grasping tips




SF-11226-1 1.8 MICS Capsulorhexis Forceps
1.8mm incision



SF-11173-1 Capsulorhexis Forceps
85mm, Straight shafts, very delicate grasping tips



SF-11173-1 1.8 MICS Capsulorhexis Forceps
1.8mm incision




SF-11173-2 Capsulorhexis Forceps
105mm, Straight shafts, very delicate grasping tips



SF-11173-2 1.8 MICS Capsulorhexis Forceps
1.8mm incision


FORCEPS



SF-11224-1 Capsulorhexis Forceps
105mm, Straight shafts, very delicate grasping tips




SF-11224-1 1.8 MICS Capsulorhexis Forceps
1.8mm incision



SF-11201-1 Capsulorhexis Forceps
105mm, Straight shafts, very delicate grasping tips




SF-11201-1 1.8 MICS Capsulorhexis Forceps
1.8mm incision





SF-11225-1 Capsulorhexis Forceps
120mm, Straight shafts, very delicate grasping tips



SF-11225-1 1.8 MICS Capsulorhexis Forceps
1.8mm incision




SF-32130-1 MICS Capsulorhexis Forceps
115mm, Cross action, curved shafts, very delicate grasping tips, 1.8mm incision




SF-11011-1 Tying Forceps
85mm, straight shafts with tying platforms,







SF-11011-2 Tying Forceps
105mm, straight shafts with tying platforms,




SF-11031-1 Tying Forceps
105mm, straight shafts with tying platforms,



SF-11018-1 Tying Forceps
85mm, curved shafts with tying platforms,




SF-11018-6 Tying Forceps
105mm, curved shafts with tying platforms,




SF-11018-4 Tying Forceps
105mm, curved shafts with tying platforms,




FORCEPS




SF-11013-1 Mcpherson Tying Forceps:
85mm, angled shafts with tying platforms,




SF-11013-2 Mcpherson Tying Forceps
105mm, angled shafts with tying platforms,




SF-11013-5 Mcpherson Tying Forceps
105mm, angled shafts with tying platforms,




SF-11282-1 Lens Forceps
115mm, curved shafts with 1.0mm ring tips,




SF-11282-2 Lens Forceps
102mm, angled shafts with 1.0mm round tips,




SF-11331-1 Lens Implantation Forceps
105mm, angled shafts with smooth jaws,




SF-11321-1 Mcpherson Lens Forceps
105mm, angled shafts with smooth jaws,



SF-11398-1 Akahoshi Prechopper:
112mm, Maximum tip opening 3.0mm




SF-31401-1 Akahoshi Prechopper:
116mm, straight shafts




Akahoshi Prechopper:
SF-11402-1 straight shafts 116mm
SF-31403-1 angled shafts 113mm




FORCEPS




SF-11400-1 Conjunctiva Forceps:
109mm, straight, 1mm ring tip




SF-15035-1 Perone LASIK Flap Forceps
109mm, curved shafts 2mm ring tip



SF-15036-1 Perone LASIK Flap Forceps
109mm, curved shafts, 1.7mm x 1.1mm tip




SF-11415-1 Chalazion Forceps
100mm, round 13mm




SF-11417-1 Chalazion Forceps
89mm, oval 24mm x 16mm







SF-11405-1 Chalazion Forceps
100mm, round 14.5mm




SF-11406-1 Chalazion Forceps
90mm, round 11mm



SF-11407-1 Chalazion Forceps
105mm, round 17mm



SF-11408-1 Chalazion Forceps
102mm, oval 16mm x 14.5mm



SF-11410-1 Chalazion Forceps
98mm, oval 13mm x 10.5mm



FORCEPS



SF-51181-1 Serrated Forceps
85mm, angled shafts with serrated tips



SF-51160-1 Lenticule Extraction Forceps
116mm, Cross action handle with serrated tips



SF-51070-1 Meibomian forceps
108mm




SF-11419-2 Compressing Forceps
115mm



SF-11414-1 Moorfields Forceps
110mm, serrated tips



SF-11412-1 Ciliary Forceps
72mm



SF-11416-1 Scleral Plug Forceps
113mm



SF-11421-1 Jewellers Forceps
115mm, straight, 0.05mm x 0.01mm tip



SF-11423-1 Jewellers Forceps
115mm, angled, 0.05mm x 0.01mm tip



SF-11424-1 Jewellers Forceps
115mm, curved, 0.07mm x 0.03mm tip

HOLDERS



SH-32008-1 Needle Holder
115mm, 9mm curved jaws, flat handle without lock



SH-12205-1 Needle Holder
140mm, 10mm straight jaws, flat handle with lock



SH-12204-1 Needle Holder
140mm, 10mm straight jaws, flat handle without lock



SH-12214-1 Needle Holder
115mm, 9mm straight jaws, flat handle with lock



SH-12213-1 Needle Holder
115mm, 9mm straight jaws, flat handle with lock



SH-12016-1 Needle Holder
115mm, 11mm curved jaws, round handle with lock



SH-12011-1 Needle Holder
115mm, 11mm curved jaws, round handle without lock



SH-12192-1 Needle Holder
140mm, 9mm curved jaws, flat handle with lock



SH-12191-1 Needle Holder
140mm, 9mm curved jaws, flat handle without lock



SH-12212-1 Needle Holder
115mm, 9mm curved jaws, flat handle with lock

HOLDERS



SH-12211-1 Needle Holder
115mm, 9mm curved jaws, flat handle with lock



SH-12243-1 Needle Holder
115mm, 9mm curved jaws, round handle with lock



SH-12235-1 Needle Holder
115mm, 9mm curved jaws, round handle without lock



SH-12177-1 Needle Holder
115mm, 9mm curved jaws, flat handle with lock



SH-12176-1 Needle Holder
115mm, 9mm curved jaws, flat handle without lock

HOLDERS

SH-35015-1 Towel Clamps
115mm, flat handle




SH-11414-1 Mosquito Forceps
125mm, straight 24mm serrated jaws,



SH-11415-1 Mosquito Forceps
125mm, curved 24mm serrated jaws,



SH-20000-1 Bulldog Clips
40mm, flat handle




Ramplays Sponge Holding Forceps:
SH-20003-1: 190mm
SH-20004-1: 145mm



SCISSORS



SH-30004-1 3# Blade Handle
125mm, from 0cm-- 10cm



SS-13080-1 Lens Cutter
85mm, angled, sharp tips, blades length 6mm



SS-13075-1 Castroviejo Scissors
115mm, curved, blunt tips, tip to pivot length 14.5mm



SS-13076-1 Castroviejo Scissors
115mm, straight, blunt tips, tip to pivot length 14.5mm



SS-13026-1 Conjunctiva Scissors
105mm, curved, blunt tips, tip to pivot length 10.5mm



SS-13071R-1 Conjunctiva Scissors
110mm, curved, blunt tips, tip to pivot length 12.5mm



SS-13056-1 Osher Universal Scissors
120mm, curved, blunt tips, tip to pivot length 16mm



SS-13030R-1 Castroviejo Corneal Scissors
115mm, curved, blunt tips, tip to pivot length 16mm



SS-13031-1 Castroviejo Corneal Scissors
115mm, curved, blunt tips, tip to pivot length 16mm



SS-13028-1 Castroviejo Corneal Scissors
115mm, curved left, blunt tips, tip to pivot length 13.5mm

SCISSORS



SS-13027-1 Castroviejo Corneal Scissors
115mm, curved right, blunt tips, tip to pivot length 13.5mm



SS-13044-1 Castroviejo Corneal Scissors
115mm, curved left, blunt tips, tip to pivot length 11.5mm



SS-13043-1 Castroviejo Corneal Scissors
105mm, curved right, blunt tips, tip to pivot length 11.5mm



SS-13073R-1 Castroviejo Corneal Scissors
110mm, curved, blunt tips, tip to pivot length 11.5mm



Castroviejo Corneal Scissors: curved, blunt tips
curved, blunt tips tip to pivot length 16mm
SS-13072R-1 115mm
SS-13060-1 135mm

SCISSORS



SS-13029-1 Conjunctiva Corneal Scissors
115mm, angled, blunt tips, blades length 10mm



SS-13042-1 Capsulotomy Vannas Scissors
117mm, curved, sharp tips, tip to pivot length 16mm



SS-13002-1 Vannas Scissors
81mm, curved, sharp tips, tip to pivot length 10mm



SS-13001-1 Vannas Scissors
81mm, straight, sharp tips, tip to pivot length 10mm



SS-13007-1 Vannas Scissors
78mm, angled, sharp tips, tip to pivot length 9mm

SCISSORS



SS-13033-1 Iris Scissors
80mm, curved right, sharp tips, cutting 4.0mm




SS-13034-1 Iris Scissors
110mm, curved right, sharp tips, blade length 7.5mm



SS-13032-1 Iris Scissors
82mm, straight, sharp tips, cutting 5.5mm



SS-13045-1 Iris Scissors
112mm, straight, sharp tips, blade length 9mm



SS-13074R-1 Westcott Stitch Scissors
112mm, curved, sharp tips, tip to pivot length 13.5mm

SCISSORS



SS-13078-1 Westcott Stitch Scissors
115mm, curved, sharp tips, tip to pivot length 16mm



SS-13059-1 Westcott Stitch Scissors
115mm, curved, sharp tips, tip to pivot length 12mm



SS-13037-1 Vannas Capsulotomy Scissors
87mm, curved, sharp tips, blades length 10.5mm



SS-13038-1 Vannas Capsulotomy Scissors
88mm, curved, sharp tips, blades length 12mm



SS-13035-1 Vannas Capsulotomy Scissors
112mm, straight, sharp tips, blades length 11mm

SCISSORS



SS-13036-1 Vannas Capsulotomy Scissors
115mm, straight, sharp tips, blades length 13mm



SS-13039-1 Vannas Capsulotomy Scissors
115mm, angled, sharp tips, blades length 11mm



SS-13040-1 Vannas Capsulotomy Scissors
115mm, angled, sharp tips, blades length 12.5mm




Tenotomy Scissors: curved, blunt tips /
SS-13079-1: 115mm, tip to pivot length 16.5mm / 115mm,
SS-13058-1: 125mm, tip to pivot length 20mm / 125mm,




SS-13065-1 Westcott Stitch Scissors
115mm, curved, blunt tips, tip to pivot length 16.5mm

SCISSORS



SS-13055-1 Westcott Stitch Scissors
110mm, curved, sharp tips, tip to pivot length 21mm



SS-13055-2 Westcott Stitch Scissors
110mm, curved, blunt tips, tip to pivot length 21mm



SS-13067-1 Vannas Scissors
81mm, straight, sharp tips, blades length 7.5mm



SS-13066-1 Vannas Capsulotomy Scissors
85mm, angled, sharp tips, blades length 9mm



SS-13062-1 General Scissors
128mm, straight, blunt tips



SCISSORS




SS-13062B-1 General Scissors
89mm, straight, sharp tips




SS-13064B-1 General Scissors
89mm, straight, blunt tips




SS-31286-1 General Scissors
98mm, straight, blunt tips



Iris Hook:
SB-14101H-1: H-shaped tip, straight / H
SB-14101H-2: H-shaped tip, angled / H




Kuglen Iris Hook:
SB-14101-1: K-shaped tip, angled / K
SB-14102-1: K-shaped tip, straight / K



Jaffe-Knolle Iris Hook:
SB-14108-1: J-shaped tip, angled / J
SB-14109-1: J-shaped tip, straight / J



SB-14143-1 Pupil Snapper Hook
With smoothly polished notch and sharp inner edge




Lens Manipulator:
SB-14104-1: T-shaped tip, straight / T
SB-14103-1: T-shaped tip, angled / T



Lens Manipulator: Y nucleus rotator, angled

SB-14105-1 Vertical
SB-14105H-1 Horizontal





Lens Manipulator:
SB-14145-1: V-shaped tip, angled / V
SB-14144-1: V-shaped tip, straight / V

PROBES, HOOKS & SPATULAS



Lens Manipulator:
SB-14106-1: Angled
SB-14107-1: Straight



SB-14110-1 Sinskey Hook
0.7mm tip, angled


SB-34077-1 Reverse Sinskey Hook
0.7mm tip, angled



SB-14111-1 Lens Manipulator
1.1mm tip, angled



SB-14113-1 Drysdale Manipulator
Paddle shaped tip, angled




SB-14135-1 Nucleus Manipulator
Angled


PROBES, HOOKS&SPATULAS




SB-14120-1 Nucleus Manipulator
Connor wand, straight



SB-14140-1 Nucleus Manipulator
Connor wand, angled



SB-14150-1 Nucleus Manipulator
Connor wand, angled, oval tip



SB-14151-1 Lens Manipulator
Connor wand, angled




SB-14112-1 Lens Manipulator
Straight



PROBES, HOOKS&SPATULAS




Lens Manipulator
SB-34125-1 Angled shaft
SB-34126-1 Straight shaft




SB-14127-1 Lens Manipulator
Minami "M" Hook



SB-14A70-1 ICL Manipulator
Discs tip



SB-1905B-1 Lens Manipulator
Angled



SB-1709A-1 Lens Manipulator
Angled with the stage at the neck

PROBES, HOOKS & SPATULAS




SB-16000-1 Cautery
5mm x 7mm oval ball with 6mm long tip



SB-14181-1 Depressor
Two 10mm disc end, 92mm



Foreign Body Remover
SB-14182-1 Straight
SB-14182-2 Curved




SB-14121-1 Kuglen Iris Hook
"S" shaped,




SB-14121-2 Kuglen Iris Hook
Curved 1.0mm hook




PROBES, HOOKS&SPATULAS




SB-34185-1 Double ended Lens Manipulator



SB-11297-1 Lenticule Separator
Pocketing hook and spatula tip curved separator



SB-14087-1 Double ended Manipulator
One side cylinder, one side Y shape manipulator




SB-14095-1 Double ended Manipulator
One side disc, one side mushroom manipulator




SB-14301-2 Double ended Kansas Nucleus Vectis


PROBES, HOOKS & SPATULAS




SB-14302-1 Double ended Chopper and Manipulator
Paddle manipulator, 60° chopper





SB-14303-1 Double ended Chopper and Manipulator
Connor manipulator, sharp chopper



SB-16002-1 Double ended Mushroom Manipulator and Phaco Chopper



SB-16003-1 Double ended Manipulator and Phaco Chopper
V shape manipulator, 90° chopper




SB-44002-1 Femto Flap Lifter & Retreatment Spatula




PROBES, HOOKS & SPATULAS




Phaco Chopper
1.2mm curved tip,
SB-14122-1 angled
SB-14321-1 straight




Phaco Chopper
1.1mm x 0.25mm tip 45°
SB-14125L-1 left hand
SB-14125-1 right hand



Nichamin Phaco Chopper
1.2mm x 0.25mm tip 90°
SB-14123-1 angled
SB-14124-1 curved




Phaco Chopper
1.6mm x 0.30mm tip 45° curved
SB-14126L-1 left hand
SB-14126-1 right hand




SB-34025-1 Nagahara Phaco Chopper
1.6mm x 0.30mm tip 60°



PROBES, HOOKS & SPATULAS




SB-14128-1 Phaco Chopper
1.6mm x 0.3mm tip 90°




SB-14129-1 Phaco Chopper
1.6mm x 0.3mm curved tip 90°



SB-14130-1 Claw Phaco Chopper
Triple inner blades



Phaco Chopper
1.6mm x 0.3mm tip 90°
SB-14136-1 angled
SB-14137-1 curved




Phaco Chopper
1.6mm x 0.3mm tip 45°
SB-14138L-1 left hand
SB-14138-1 right hand




PROBES, HOOKS & SPATULAS



SB-14139-1 Haefliger Phaco Chopper
Chopper and manipulator in one




SB-14141-2 Phaco Chopper
Ambidextrous Chopper



SB-14416-1 Phaco Chopper
Ball tip




SB-14131-1 Nealiann Nucleus Divider
Ambidextrous Chopper



SB-15027-1 Corneal Epithelium Spatula




PROBES, HOOKS&SPATULAS




SB-34022-2 Rosen Chopper
1.0mm x 0.3mm straight tip




SB-14116-1 Kansas Nucleus Bisector
For manual small incision surgery



SB-34092-1 Phaco Chopper
Angled




SB-14168-1 Epi-peeler




SB-14169-1 Hoe Type Epi-peeler



PROBES, HOOKS&SPATULAS



SB-14161-1 Helveston Scleral Marking Ruler
0~10mm



Trabeculectomy Punch
SB-20011-1 0.75mm
SB-20011-2 1.00mm
SB-20011-3 1.50mm
SB-20011-4 2.00mm



SB-35074-1 Scleral Depressor
Single bulbous end



SB-35075-1 Scleral Depressor
Single hammer end



SB-35076-1 Scleral Depressor
Double ended




PROBES, HOOKS & SPATULAS




SB-14118-1 Lens Loop
5mm x 7mm



SB-14118-2 Lens Loop
3mm x 8mm



Repositor
SB-14114-1 0.7mm
SB-14115-1 0.5mm



SB-34110-1 Double ended Repositor



Curette
SB-34008-1 1.0mm
SB-34006-1 1.5mm
SB-34007-1 2.5mm
SB-34004-1 3.0mm
SB-34009-1 3.5mm



PROBES, HOOKS & SPATULAS



Curette
SB-34055-1 Oval shape



Double ended Curette
SB-34003-1 1.0mm / 2.0mm
SB-34179-1 1.5mm / 2.0mm
SB-34127-1 2.0mm / 3.0mm



SB-34186-1 Double ended Curette
One side 2.0mm round tip, the other side oval tip 2.0mm




SB-15032-2 Fine Thornton Fixation Ring



Fine Thornton Fixation Ring
SB-14148-1 Swivel head with 16mm Dia. 14mm cut out
SB-14148-2 Swivel head with 13mm Dia. 11mm cut out
SB-14149-1 Swivel head with 12mm Dia. 9mm cut out

PROBES, HOOKS & SPATULAS



Lasek Alcohol Well

SB-14171-1 6.0mm

SB-14173-1 8.0mm

SB-14171-2 6.5mm


SB-14173-2 8.5mm

SB-14172-1 7.0mm

SB-14174-1 9.0mm


SB-14172-2 7.5mm

SB-14174-2 9.5mm




SB-14133A-1 Retinal Hook

0.5mm tip




SB-14133B-1 Retinal Hook

Round tip



SB-14133C-1 Retinal Hook


110° angled tip




SB-14133D-1 Retinal Hook




PROBES, HOOKS & SPATULAS




SB-14176-1 Intra-Lasik Flap Separator



SB-14119-1 Retinal Detachment Hook
1.5mm oval hole



Jameson Strabismus Hook
SB-14170-1 6.5mm
SB-14170-2 8mm



SB-14175-1 Chavasse Hook
12mm



SB-14134-1 Freer Lacrimal Chisel
5.5mm x 12mm tip

PROBES, HOOKS & SPATULAS



SB-14166-1 Arruga Globe Retractor
6.6mm wide & 6.6mm wide



SB-14167-1 Arruga Globe Retractor
10.5mm wide & 3.0mm wide

SPECULUMS

Jaffe Lid Retractors x 2

对拉开睑器 (2只)


SO-15077-1 16.5mm blades, overall length 55mm

SO-15077-2 12mm blades, overall length 55mm




SO-15001-1 Wire Speculum

15.5mm blades, overall length 48mm




SO-15007-1 Kratz Barraquer Wire Speculum

15.5mm blades, overall length 48mm




SO-15009-1 Naqahara Wire Speculum

15.5mm blades, overall length 48mm



SO-15010-1 Grandon Wire Speculum

13mm blades, overall length 46mm




SPECULUMS


SO-15013-1 Wire Speculum
15.5mm blades, overall length 50mm




SO-15019-1 Wire Speculum
15mm blades, overall length 48mm



SO-15020-1 Wire Speculum
16mm blades, overall length 48mm




SO-15021-1 Pierce Type Adjustable Speculum
16mm blades, overall length 30mm




Wire Speculum

SO-15017-1 12mm blades, overall length 46mm
SO-15018-1 16mm blades, overall length 46mm




SPECULUMS

SO-15051-1 Wire Speculum
16mm blades, overall length 40mm



Wire Speculum


SO-15053-1 11mm blades, overall length 40mm
SO-15052-1 11mm blades, overall length 42mm



SO-15055-1 Wire Speculum
16mm blades, overall length 44mm



SO-15056-1 Wire Speculum
11mm blades, overall length 42mm




SO-15090N-1 Wire Speculum
13mm blades, overall length 42mm



Wire Speculum

SO-15059-1 10.5mm blades, overall length 30mm

SO-15060-1 9mm blades, overall length 29mm



SO-15075-1 Wire Speculum with Discs
8.5mm and 10.5mm blades, overall length 50mm



Sprung Speculum

SO-15078-1 6mm Paediatric

SO-15078-2 4mm Neonate



SO-15015-1 Barraquer Solid Speculum
14mm blades, overall length 42mm




SO-15015-2 Barraquer Solid Speculum
13.5mm blades, overall length 45mm



SPECULUMS

SO-15008-1 Barraquer Solid Speculum
14.5mm blades, overall length 45mm



SO-35002-1
Sliding Speculum



SO-15002-1 Lieberman Adjustable Speculum
V-shaped Wire, overall length 80mm



SO-15002N-1 Adjustable Speculum
15mm blades, overall length 80mm



SO-15005-1 Adjustable Speculum
15mm blades, overall length 80mm




SPECULUMS

SO-15005N-1 Reversible Adjustable Speculum
15mm blades, overall length 80mm




SO-15014-1 Adjustable Speculum
Swivel blades, overall length 80mm



SO-31431-1 Adjustable Speculum
19mm blades, overall length 80mm



SO-15063-1 Adjustable Speculum
15mm blades, overall length 80mm



SO-15006-1 Adjustable Speculum
12.5mm blades, overall length 78mm



SPECULUMS

SO-15011-1 Adjustable Speculum
13mm blades, overall length 76mm



SO-15022-1 Adjustable Speculum
solid blades, overall length 78mm



SO-15029-1 Adjustable Speculum
13mm blades, overall length 78mm



SO-15033-1 Adjustable Speculum
14mm blades, overall length 78mm




SO-15034-1 Adjustable Speculum
14mm blades, overall length 78mm



SPECULUMS

SO-15004-1
Lindstrom-chu Aspirating Adjustable Speculum
overall length 78mm



SM-34143-1 Braunstein Fixed Caliper
3.5mm / 4.0mm




SM-34143-2 Braunstein Fixed Caliper
3.5mm / 4.0mm




SM-20005-1 Castroviejo Caliper
Calibrated every 0.5mm from 0 to 20mm




MARKERS



SM-34001-1 Mendez Degree Gauge
Calibrated every 10° from 0° to 180°
External dia 19.5mm, internal dia 12.5mm



SM-34002-1 Optical Zone Marker
Target center, 4.0mm dia.



SM-34002-2 Optical Zone Marker
Target center, External dia 10mm, internal dia 7mm




SM-34004-1 LASIK Corneal Marker
Two blunt fins




SM-34003-1 Scleral Marker
3.5mm / 4.0mm




Vitreoretinal Scissors




SV-13151-1 20G
 SV-13151-2 23G
 SV-13151-3 25G
 Straight Scissors, sharp tips,
 straight tube




SV-13152-1 20G
 SV-13152-2 23G
 SV-13152-3 25G
 Straight Scissors, blunt tips,
 straight tube




SV-13153-1 20G
 SV-13153-2 23G
 SV-13153-3 25G
 Curved Scissors, sharp tips
 straight tube



SV-13154-1 20G
 SV-13154-2 23G
 SV-13154-3 25G
 Angled Scissors, sharp tips
 straight tube



SV-13155-1 20G
 SV-13155-2 23G
 SV-13155-3 25G
 Vertical Scissors, 40°, straight tube



SV-13156-1 20G
 SV-13156-2 23G
 SV-13156-3 25G
 Vertical Scissors, 80°, straight tube



SV-13169-1 20G
 SV-13169-2 23G
 SV-13169-3 25G
 Straight Scissors, sharp tips
 28mm curved tube




SV-13171-1 19G
 IOL Scissors, straight tube

SV-13170-1 20G
 SV-13170-2 23G
 SV-13170-3 25G
 Straight Scissors, sharp tips
 15mm curved tube


VITREORETINAL



Vitreoretinal Forceps




SV-11502-1 20G
Cap Jaws Forceps, straight tube




SV-11503-1 20G
SV-11503-2 23G
SV-11503-3 25G
Capsulorhexis Forceps, curved tube




SV-11504-1 20G
SV-11504-2 23G
SV-11504-3 25G
Capsulorhexis Forceps, straight tube




SV-11505-1 20G
SV-11505-2 23G
SV-11505-3 25G
Gripping Forceps, sharp tips,
straight tube




SV-11506-1 20G
SV-11506-2 23G
SV-11506-3 25G
Gripping Forceps, round tips,
straight tube




SV-11507-1 20G
SV-11507-2 23G
SV-11507-3 25G
Serrated Forceps, straight tube



SV-11507C-1 20G
SV-11507C-2 23G
SV-11507C-3 25G
Serrated Forceps, curved tube



SV-11508-1 20G
SV-11508-2 23G
SV-11508-3 25G
Toothed Forceps, 1 x 2 teeth,
straight tube



SV-11509-1 20G
SV-11509-2 23G
SV-11509-3 25G
End Gripping Forceps, straight tube


VITREORETINAL



Vitreoretinal Forceps



SV-11510-1 20G
SV-11510-2 23G
SV-11510-3 25G
Angled Forceps, smooth, straight tube



SV-11511-1 20G
SV-11511-2 23G
SV-11511-3 25G
Angled Forceps, grooved, straight tube



SV-11512-1 20G
SV-11512-2 23G
SV-11512-3 25G
End Grasping Forceps, straight tube



SV-11520-1 20G
SV-11520-2 23G
SV-11520-3 25G
Asymmetrical forceps, straight tube




SV-11578-1 20G
SV-11578-2 23G
ICL forceps Ribbed, straight tube




SV-11580-1 20G
ICL forceps, straight tube



SV-11581-1 20G
SV-11581-2 23G
SV-11581-3 25G
Eckardt End Gripping Forceps, straight tube



SV-30001-1 20G
Foreign Body Forceps
2 jaws



SV-11528-1 20G
Foreign Body Forceps
3 jaws


CANNULAS



SC-34001-1 Infusion Handpiece
7mm handle, overall length 90mm



SC-34002-1 Infusion Handpiece
7mm handle, double ends, overall length 95mm



SC-34003-1 Lacrimal Probe
4pcs, overall length 142mm



Pigtail Lacrimal Probe
SC-34004-1 6mm
SC-34004-2 7mm
SC-34004-3 8mm
SC-34004-4 9mm



SC-34081-1 Punctal Dilator
sharp 17mm taper

CANNULAS



SC-34013-1 Punctal Dilator
sharp 15mm taper



SC-34015-1 Punctal Dilator
double ends,
sharp 17mm & 13mm taper



SC-34021-1 LASIK Cannula
Angled shaft 17mm from tip to bend, 25g



SC-34023-1 Posterior Capsule Polisher
10mm flatted angled and sandblasted tip with
irrigation holes on tip, 18g



SC-34023-2 Posterior Capsule Polisher
3mm flatted angled and sandblasted tip with
irrigation holes on tip, 21g

CANNULAS



SC-34023-3 Posterior Capsule Polisher
1.5mm smooth bulbous and sandblasted tip with irrigation holes on tip, 23g

Irrigating Vecits 18g

- SC-34024-1 4mm x 6mm loop
- SC-34024-2 3mm x 8mm loop
- SC-34024-3 5mm x 7mm loop
- SC-34024-4 3mm x 12mm loop

4mm x 6mm



3mm x 8mm



5mm x 7mm



3mm x 12mm



Irrigating- Aspirating Cannula

- SC-34025-1 23g 6
- SC-34025-2 22g 7
- SC-34025-3 21g 8




Irrigating- Aspirating Cannula

- SC-34025-4 22g, angled, left 7
- SC-34025-5 22g, angled, right 7



- Lewicky Anterior Chamber Maintainer
18g, 30mm silicone tube
SC-34028-1



CANNULAS

SC-10026-1 Lacrimal Cannula
1 x 5, 25g, 23g, 22g, 21g, 20g



SC-34029-1 Lacrimal Cannula
45° angled 10mm tip, 23g



SC-34029-2 Lacrimal Cannula
straight, 25g



SC-34029-3 Lacrimal Cannula
curved tip, 23g




SC-10028-1 LASIK Cannula
plastic, single hole at tip




PACK

PST 01 Sub-Tenons Pack
SO-15019-1
SS-13065-1
SF-11414-1



PST 02 Sub-Tenons Pack
SO-15019-1
SS-13065-1
SF-11061-2



PST 03 Suture Pack
SS-13062B-1
SH-32008-1
SF-11011-2



PST 04 Suture Pack
SS-13062B-1
SH-32008-1
SF-11018-1



PST 05 Cataract Set
SO-15005N-1
SF-11201C-1
SF-11111-2
SB-14078-1



PST 06 Phaco Set
SO-15005N-1
SF-11201-1
SS-31286-1
SB-14095-1



